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INFORMATION

Date of last update: 2025-03-12

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

This document applies to the following products: F21 Maverick Xtreme AT, F22 Maverick Comfort AT

The prosthetic foot has fibreglass spring elements. The heel stiffness can be increased with heel wedges.

2 Safety-related information

CAUTION! Risk of injury and risk of product damage

- ▶ Observe the information on product combinations in the instructions for use for the products.
- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 15).
- ▶ To avoid the risk of injury and product damage, do not use the product beyond the tested lifetime.
- ▶ To avoid the risk of injury and product damage, only use the product for a single patient.
- ▶ To prevent mechanical damage, treat the product with care.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

NOTICE! Danger of product damage and restricted functionality

- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Clean the product if it gets dirty. If necessary, have it inspected, repaired or replaced by the manufacturer or a specialist workshop.

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover behaviour) or delamination of the spring are indications of loss of functionality. Unusual noises can indicate a loss of functionality.

3 Intended use

3.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

3.2 Operating conditions

The MOBIS classification describes the mobility grade and body weight, and makes it easy to identify compatible components.

Maverick Xtreme AT

The product is recommended for mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands).

Maverick Comfort AT

The product is recommended for mobility grade 2 (restricted outdoor walker) and mobility grade 3 (unrestricted outdoor walker).

The table that follows shows the suitable spring stiffness of the prosthetic foot, matching the patient's body weight and activity.

Spring stiffness relative to body weight and activity level			
Body weight [kg]	Low	Normal	High
44 to 52	1	1	2
53 to 59	1	2	3
60 to 68	2	3	4
69 to 77	3	4	5
78 to 88	4	5	6
89 to 100	5	6	7
101 to 116	6	7	8
117 to 130	7	8	9
131 to 147	8	9	–
148 to 166	9	–	–

3.3 Environmental conditions

Storage and transport

Storage temperature: -20 °C to +60 °C, relative humidity: 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions

Temperature during use: -10 °C to +45 °C

Chemicals/liquids: Fresh water, salt water, perspiration, urine, soapsuds, chlorine water

Moisture: Submersion: max. 1 h in 3 m depth, relative humidity: no restrictions

Solids: Dust, occasional contact with sand

Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and tear and damage (see page 20).

Prohibited environmental conditions
Solids: Highly liquid-binding particles (e.g. talcum), continuous contact with sand
Chemicals/liquids: Acids, continuous use in liquid media

3.4 Product combinations

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Limited combination options for Ottobock components

The prosthetic foot generates high moments in the ankle area. Use structural components with higher weight limits:

Body weight (high activity) [kg]	up to 59	up to 77		up to 100		up to 130
Foot size [cm]	up to 26	up to 28	from 29	up to 28	from 29	up to 31
Structural component weight limit [kg]	≥75	≥100	≥125	≥125	≥150	≥150

Body weight (normal activity) [kg]				up to 100		up to 130
Foot size [cm]				from 29		up to 31
Structural component weight limit [kg]				≥125		≥150

3.5 Lifetime

Prosthetic foot

Depending on the user's activity level, the maximum lifetime of the product is 3 years.

Footshell, protective sock

The product is a wear part, which means it is subject to normal wear and tear.

4 General safety instructions

5 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	–
1	Prosthetic foot	–
1	Protective sock (black)	SL=Spectra-Sock-7
1	Heel wedge set	2F11=1

Spare parts/accessories	
Designation	Reference number
Footshell (without cap)	FTC-2A-1*
Tool to change the footshell	ACC-00-10300-00

6 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- Observe the alignment, assembly and adjustment instructions.

NOTICE

Altering the prosthetic foot or footshell

Premature wear resulting from damage to the product

- Do not alter the prosthetic foot or footshell.

INFORMATION

A plastic part protects the connecting section of the product from scratches during the alignment and fitting.

- Remove the plastic part before the patient leaves the fitting area.

6.1 Applying/removing the footshell

- Pull the protective sock over the prosthetic foot to prevent noises in the footshell.
- **CAUTION! Always use the prosthetic foot with the footshell.**
Apply or remove the footshell as described in the footshell instructions for use.

6.2 Bench Alignment

Bench alignment, transtibial prosthesis

Bench alignment process

Required materials: 662M4 goniometer, 743S12 heel height measuring device, 743A80 50:50 gauge, alignment tool (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly)

Assemble and align the prosthetic components in the alignment apparatus according to the following specifications:

1

Bench alignment process	
	Dorsiflexion/plantar flexion: Adjust the prosthetic foot to the heel height of the shoe. (The prosthetic foot is designed for shoes with a heel height of 10 mm .)
②	Adduction/abduction: Adjust the angle of the prosthetic socket in the frontal plane.
③	Flexion/extension: Adjust the angle of the prosthetic socket in the sagittal plane.
④	Linear shifting of the prosthetic socket: Shift the prosthetic socket so that the plumb line runs along the front edge of the adapter on the prosthetic foot (see fig. 3).

Bench alignment, thigh/knee disarticulation prosthesis

- Observe the instructions for use for the prosthetic knee joint.

6.3 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.

6.4 Dynamic Trial Fitting

The optimal gait pattern is developed during dynamic fitting. The prosthesis is aligned optimally in the frontal plane and the sagittal plane for this purpose.

- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal planes is achieved after the leg begins to bear weight following the heel strike. Avoid medial movement of the knee joint.
 - If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction.
 - If the knee joint moves in the medial direction in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.
- After completing the dynamic fitting and the walking exercises: Remove the plastic part from the pyramid.

Heel too soft	
Symptoms	Possible solutions
• Foot makes full ground contact too quickly • Forefoot feels too stiff • Knee goes into hyperextension	• Shift prosthetic socket forward relative to the foot • Use a heel wedge

Heel too hard	
Symptoms	Possible solutions
• Fast knee flexion, low stability	• Shift prosthetic socket back relative to the foot

Heel too hard	
• Transition from heel strike to toe-off too fast • Energy return feels low	• Reduce the heel stiffness (shift or remove heel wedge)

Prosthetic foot too stiff	
Symptoms	Possible solutions
• Minimal rollover movement of the prosthetic foot at low walking speed (long full ground contact)	• Select prosthetic foot in lower stiffness

Prosthetic foot too soft	
Symptoms	Possible solutions
• Clicking noise at the start of ground contact. • Very pronounced deformation of the forefoot at high activity levels	• Select prosthetic foot in higher stiffness

6.4.1 Optimising the heel characteristics

The behaviour of the prosthetic foot at heel strike and during heel contact in the mid-stance phase can be adapted with the placement of a heel wedge. The heel wedge is attached with adhesive tape for testing. For definitive assembly, it is glued to the prosthetic foot.

Fitting

- 1) Use the adhesive tape on the underside of the heel wedge.
- 2) Place the heel wedge on the base spring in the recommended position (4602824).
- 3) Set the desired stiffness by shifting the heel wedge (anterior = harder, posterior = softer).

Definitive assembly

- 1) Remove the adhesive tape from the heel wedge. This can be done using a cloth with acetone.
- 2) Roughen the contact surface on the prosthetic foot slightly using sandpaper. Remove the dust after sanding.
- 3) Glue the heel wedge to the prosthetic foot with cyanoacrylate adhesive.

7 Cleaning

> **Allowable cleaning agent:** pH-neutral soap (e.g. 453H10 Derma Clean)

- 1) **NOTICE! To avoid product damage, only use the allowable cleaning agents.**

Clean the product with clear fresh water and a pH neutral soap.

- 2) Remove debris from water drainage contours using a toothpick and rinse.
- 3) Rinse the soap away with clear fresh water. In doing so, rinse the foot-shell until all debris has been removed.
- 4) Dry the product with a soft cloth.
- 5) Allow to air dry in order to remove residual moisture.

8 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ During regular inspection: Inspect the prosthesis for signs of wear and check its function.
- ▶ Inspect the product for visible damage every six months. Make additional inspection appointments as needed (e.g. for highly active or heavy users).

9 Disposal

Do not dispose of the product in household waste. Improper disposal is harmful to health and the environment. Observe the instructions of the local authority regarding return, collection and disposal.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

11 Technical data

F21 Maverick Xtreme AT									
Sizes [cm]	23	24	25	26	27	28	29	30	
Heel height [mm]	10								

F21 Maverick Xtreme AT								
Sizes [cm]	23	24	25	26	27	28	29	30
System height [mm]	144			152			165	
Build height [mm]	162			170			183	
Average weight with footshell [g]	665			820			910	
Max. body weight [kg]	116	130		147	166			
Mobility grade	3, 4							

F22 Maverick Comfort AT								
Sizes [cm]	23	24	25	26	27	28	29	30
Heel height [mm]	10							
System height [mm]	98			99			101	
Build height [mm]	116			117			119	
Average weight with footshell [g]	615			750			830	
Max. body weight [kg]	116	130		147	166			
Mobility grade	2, 3							

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