12K41, 12K42, 12K44, 12K50

DE Gebrauchsanweisung (Benutzer)	4
EN Instructions for use (user)	14
FR Instructions d'utilisation (utilisateur)	25
IT Istruzioni per l'uso (utilizzatore)	36
ES Instrucciones de uso (usuario)	47
PT Instruções de utilização (utilizador)	58
NL Gebruiksaanwijzing (gebruiker)	69
SV Bruksanvisning (användare)	80
DA Brugsanvisning (bruger)	90
NO Bruksanvisning (bruker)	101
FI Käyttöohje (käyttäjä)	111
PL Instrukcja używania (użytkownik)	122
HU Használati útmutató (használó)	133
CS Návod k použití (uživatel)	143
TR Kullanım kılavuzu (kullanıcı)	153
EL Οδηγίες χρήσης (Χρήστης)	164
RU Руководство по применению (для пользователей)	175
UK Інструкція з використання (для користувачів)	187
JA 取扱説明書(装着者用)	198

**Order your
free printed copy**

order-ifu@ottobock.com

Document: 647H919 Version: 05

<https://product-documents.ottobock.com/IFU/INT/12K42/647H919/05/O/S/F>

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INFORMATION

Date of last update: 2024-09-16

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.


The product “Ergo Arm 12K41=* , ErgoArm plus 12K42=* , ErgoArm Hybrid plus 12K44=* , ErgoArm Electronic plus 12K50=*” is referred to as the product/elbow joint below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Design

1. Easy Plug
2. Band clamp
3. Lifter band
4. Switch cable (mechanical unlocking mechanism)
5. Handwheel
6. Humeral rotation feature adjusting screw

2.2 Function

The product is a mechanically/electronically driven elbow joint. It helps the user perform everyday tasks in combination with other prosthetic components (Combination possibilities).

Various control programmes make optimal adaptation to the individual needs and abilities possible.

The ErgoArm models differ in the following characteristics:

Models	AFB	Easy Plug	Lock	Slip-stop
ErgoArm 12K41=*	-	-	Mechanical	Mechanical
ErgoArm plus 12K42=*	x	-	Mechanical	Mechanical
ErgoArm Hybrid plus 12K44=*	x	x	Mechanical	Mechanical
ErgoArm Electronic plus 12K50=*	x	x	Mechanical or Electronic	Mechanical or Electronic

2.2.1 Definition of terms

AFB (Automatic Forearm Balance)

AFB is a mechanical flexion aid. It stores energy when the arm is extended and uses it to minimise the force required for flexion. The level of support can be adjusted to the weight of the prosthetic forearm and clothing with a handwheel.

Lock

The product is either locked mechanically or electronically (see page 15). The product can be locked and unlocked in any position, even under load, by activating the switch cable. The product has a load capacity of up to 177 N in the locked state. The lock will slip under higher loads.

Mechanical lock

The product can be manually unlocked and locked by pulling hard on the switch cable.

This function is also possible when the product is switched off or the rechargeable battery is drained. In this way, the forearm can be brought into the desired position. Mechanical unlocking is also possible under load.

Electronic lock

The lock can be released and engaged using electrode signals or switch signals (depending on the chosen switching version).

Slip-stop function

The slip-stop function makes it possible to lower the forearm in a controlled manner without having to fully release the lock.

Slip-stop function, mechanical

By pulling slightly on the switch cable (no audible click), the product unlocks so the forearm can be moved continuously. When the switch cable is released, the product is locked again.

Slip-stop function, electronic

The slip-stop function is activated by the electrode signals (depending on the chosen control version). When the electrode signal stops, the product is locked again.

3 Intended use

3.1 Indications for use

The product is to be used solely for upper limb exoprosthetic fittings.

3.2 Conditions of use

The product is suitable for users with unilateral or bilateral upper arm amputations.

The product was developed for everyday use and must not be used for unusual activities. These activities include, for example, sports/extreme sports with high strain on the wrist and/or shocks such as push-ups, downhill racing, mountain biking, free climbing, paragliding, etc.

The product is intended **exclusively** for use on one user. Use of the product by a further person is prohibited by the manufacturer.

3.3 Contraindications

- All conditions which contradict or go beyond the specifications listed in the sections on "Safety" (see page 16) and "Intended use" (see page 16).

3.4 Qualification

Treatment may only be carried out by trained, qualified personnel.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 General safety instructions

 WARNING
Operating the product near active implanted systems
Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the product.
► When operating the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
► Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

WARNING

Use of the product while operating a vehicle or machinery

- > Risk of accident due to unexpected behaviour of the product.
- > Risk of injury due to faulty control or malfunction of the product.
- ▶ Applicable national legal regulations must be complied with when operating vehicles of any kind or when operating machinery with the product.
- ▶ Have your driving ability examined and verified by an authorised test centre.
- ▶ Switch off the product before driving.

CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

CAUTION

Penetration of dirt and moisture into the product

Risk of injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.
- ▶ Do not expose the product, especially the elbow joint, to dripping or splashing water.
- ▶ Wear the product, and especially the elbow joint, under suitably resistant clothing in the rain.

CAUTION

Incorrect donning and doffing of the prosthesis

Risk of injury due to sudden flexion of the prosthesis.

- ▶ Switch off the prosthesis before donning or doffing.
- ▶ The prosthesis must always be flexed for donning and doffing.

CAUTION

Risk of pinching in the joint flexion area

Injuries due to pinching of body parts.

- ▶ Ensure that fingers and other body parts are not in this area when bending the joint.

CAUTION

Manual unlocking of elbow lock under load

Injury by release of elbow lock under load.

- ▶ Particular caution should be exercised when unlocking the elbow lock while lifting heavy loads.
- ▶ Be careful when unlocking the lock under such conditions due to the possibility of injury.

CAUTION

Mechanical stress on the product

Injury due to faulty control or malfunction of the product components.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

CAUTION

Sudden flexion of the product with high flexion assist

Risk of injury due to sudden flexion of the prosthesis.

- ▶ Observe flexion assist when taking off heavy clothing.
- ▶ Observe the flexion assist setting when the prosthesis is not used.
- ▶ The prosthesis must always be flexed for donning and doffing.

CAUTION

Removing the switch cable

Injury due to faulty control or malfunction of the product.

- ▶ The switch cable is to be used for emergency operation only. For safety reasons, it must not be removed!

CAUTION

Overheating of the product

Injury due to loss of product functionality.

- ▶ Do not expose the product to direct sunlight or heat sources.
- ▶ Avoid remaining in areas outside the allowable temperature range, specified in the section "Technical data" (Technical data).

CAUTION

Overloading due to unusual activities

Risk of injury due to unexpected product behaviour as a result of a malfunction.

- ▶ Do not expose the product to excessive strain (see page 16) or to prohibited environmental conditions (see page 24).
- ▶ Do not use the product if it is damaged or in a questionable condition (e.g. after falling). Have the product repaired or replaced by the manufacturer or a specialist workshop.

CAUTION

Signs of wear and tear on the product components

Injury due to faulty control or malfunction of the product

- ▶ In the interest of preventing injuries and to maintain product quality, we recommend performing regular service.
- ▶ Please contact your O&P professional for more information.

CAUTION

Changing prosthetic components when switched on

Risk of injury due to faulty control or malfunction of the prosthesis solution.

- ▶ Remove the rechargeable battery from the battery receptacle or switch off the prosthesis solution before changing prosthetic components (such as the terminal device).
- ▶ The 12K50=* ErgoArm can be locked and unlocked by means of the mechanical pull cable even when it is switched off.

CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ It is therefore recommended to maintain the following minimum distances from these HF communication devices:
 - Mobile phone GSM 850/GSM 900: 0.99 m
 - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
 - DECT cordless phones incl. base station: 0.35 m
 - WiFi (routers, access points ...): 0.22 m
 - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

NOTICE

Coating, gluing or painting the prosthesis

Damage or fracture due to chemical processes.

- ▶ The prosthesis must not be coated, glued or painted.

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product only as described in the section "Cleaning and care" (Cleaning and Care).

5 Scope of delivery

12K41 ErgoArm scope of delivery

- 1 pc. 12K41=* ErgoArm
- 1 pc. instructions for use (user)

12K42 ErgoArm plus scope of delivery

- 1 pc. 12K42=* ErgoArm plus
- 1 pc. instructions for use (user)

12K44 ErgoArm Hybrid plus scope of delivery

- 1 pc. 12K44=* ErgoArm Hybrid plus
- 1 pc. instructions for use (user)

12K50 ErgoArm Electronic plus scope of delivery

- 1 pc. 12K50=* ErgoArm Electronic plus
- 1 pc. instructions for use (user)



6 Use

6.1 Adjusting the flexion assist

Complete the following steps to set the flexion assist:

- 1) Flex the arm.
- 2) Use the handwheel to set the desired flexion assist.

INFORMATION: The support level can be adjusted to the weight of the clothing and the hand weight.

Position 1

Low support to allow the arm to swing freely while walking.

Position 2

Support increases progressively when flexing the arm and decreases during extension.

Position 3

Support remains constant. If set correctly, the weight of the forearm is balanced by the flexion assist and the forearm "floats".

Position 4

Low support before flexion stop.

6.2 Adjusting the humeral rotation feature

The humeral rotation feature has a bilateral stop of $\pm 80^\circ$. The resistance of the humeral rotation feature can be adjusted using the external adjusting screw.

- 1) Increase the resistance by turning it clockwise.
- 2) Reduce the resistance by turning it counter-clockwise.

7 Charging the battery

⚠ WARNING

Charging the prosthesis without taking it off

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- For safety reasons, remove the prosthesis prior to charging it.

⚠ CAUTION

Charging the product with damaged power supply unit/charger/charger cable

Falling due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

INFORMATION

Please observe the corresponding instructions for use of the charging equipment.

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product every day when used on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.

7.1 757L20 battery charger and 757B2* rechargeable battery

- 1) Plug the charging station with power supply unit into the wall socket.
→ The green LED on the power supply unit and the green LED in the middle of the LED row light up.
- 2) Insert one or two rechargeable batteries into the charging station.
- 3) The green LED for the respective charging socket lights up ● and the charging process starts.
- 4) When the rechargeable battery is fully charged, the green LED for the respective charging socket flashes ●.
- 5) Remove the rechargeable battery after the charging process is complete.
- 6) Unplug the power supply unit from the outlet.

7.1.1 Display of the current charge level

When inserting/engaging the rechargeable battery in the prosthesis, a charge level indicator is activated on the battery for a few seconds.

LED indicator	Event
●	Charge level above 50 % (illuminated in green)
● and ●	Charge level below 50 % (flashing alternating green and orange)
●	Charge level below 5 % (illuminated in orange)

7.1.2 Troubleshooting

Charging shaft	LED for the charging socket	Event
Empty	● or ● or ●	The battery charger is defective The battery charger and power supply unit must be inspected by an authorised Ottobock Service Centre.
Empty or rechargeable battery inserted	●	
Rechargeable battery inserted	●	The rechargeable battery is charging

Charging shaft	LED for the charging socket	Event
Rechargeable battery inserted		The rechargeable battery is fully charged
		The temperature of the rechargeable battery is too high Remove the rechargeable battery and let it cool down
		The inserted rechargeable battery is defective The rechargeable battery must be inspected by an authorised Ottobock Service Centre.

7.2 757L35 battery charger and 757B35=* rechargeable battery

7.2.1 Connect battery charger to the product

INFORMATION

The product cannot be used as long as the charging plug is connected to the charging receptacle. The product is deactivated for the duration of charging.

- 1) Connect the charging plug to the charging receptacle on the prosthetic socket.
→ A correct connection between the battery charger and the prosthesis is indicated by feedback (Status signals).
- 2) The charging process starts.
→ The prosthesis is switched off automatically.
- 3) Disconnect the prosthesis when the charging process is complete.

7.2.2 Safety Shutoff

The purpose of the prosthesis safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature
- Overvoltage or undervoltage
- Short circuit

After a short circuit, the charging plug must be plugged into the charging receptacle and then unplugged in order for the electronics to be activated.

7.2.3 Display of the current charge level

The charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button for less than one second.
- 2) The LED indicator on the charging receptacle provides information on the current battery charge level.

Charging receptacle	Event
	Battery fully charged (illuminated in green)
	Battery 50 % charged (illuminated in yellow)
	Battery drained (illuminated in orange)

7.2.4 Beep signals

Beep signal	Additional display	Event
1 x long	-	• Turn off prosthesis on charging receptacle • Charging begins (charging plug connected to charging receptacle) • Charging ends (charging plug disconnected from charging receptacle)
2 x short	LED on charging receptacle lights up briefly	Switch on prosthesis
3 x short	-	Battery voltage too low, prosthesis shuts off automatically

8 Cleaning

- 1) Clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock Derma Clean) when needed.
Make sure that no liquids get into the product and product components.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

9 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

The product is required for maintenance 4 weeks before, up to a maximum of 3 months after the due date.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

Components required for maintenance or repair:

The prosthesis, battery charger and power supply unit.

10 Legal information

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

10.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address:
<http://www.ottobock.com/conformity>

11 Technical data

Environmental conditions				
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85 % relative humidity, non-condensing			
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90 % relative humidity, non-condensing			
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95 % relative humidity, non-condensing			
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85 % relative humidity, non-condensing			

Reference number	12K41=*	12K42=*	12K44=*	12K50=*
Current draw (in mA)	-	-	-	3.5
Power supply	-	-	-	757B20 757B21 757B35=*
Weight (g)	%Approx. 539	%Approx. 672	%Approx. 695	%Approx. 721
Length (mm)	%Approx. 350			
Max. load in locked state (N)	177			
Flexion angle (in °)	15° – 145°			
Product behaviour during the charging process	The product is non-functional			
Expected lifetime given compliance with the recommended maintenance intervals	6 years			

12 Symbols Used

Declaration of conformity according to the applicable European directives

Serial number (YYYY WW NNN)

YYYY – year of manufacture

WW – week of manufacture

NNN – sequential number

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

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This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
T +43-1 523 37 86 · F +43-1 523 22 64
info.austria@ottobock.com · www.ottobock.com