1E57 Lo Rider

DE	Gebruuchsanweisung	3
EN	Instructions for use	8
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	19
ES	Instrucciones de uso	25
PT	Instruções de utilização	31
NL	Gebruiksaanwijzing	36
SV	Bruksanvisning	42
DA	Brugsanvisning	47
NO	Bruksanvisning	52
FI	Käyttöohje	57
PL	Instrukcja użytkowania	63
HU	Használati útmutató	68
CS	Návod k použití	74
RO	Instrucţiuni de utilizare	79
HR	Upute za uporabu	85
SL	Navodila za uporabo	90
SK	Návod na použitie	95
BG	Инструкция за употреба	101
TR	Kullanım kılavuzu	106
EL	Οδηγίες χρήσης	112
RU	Руководство по применению	118
UK	Інструкція щодо використання	124
JA	取扱説明書	129
ZH	使用说明书	134
KO	사용 설명서	139

INFORMATION

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INFORMATION

Date of last update: 2025-05-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The 1E57 Lo Rider prosthetic foot is suitable for walking on various surfaces and for a wide range of speeds.

The special combination of carbon and polyurethane absorbs heel strike. The prosthetic foot permits perceptible plantar flexion at heel strike and a natural rollover movement. Stored energy is released by the spring elements.

2 Safety-related information

CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, treat the product with care.
- ▶ Observe the information on product combinations in the instructions for use for the products.
- ▶ Observe the maximum lifetime of the product (see page 10).
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Observe the operating conditions of the product (see page 9) and do not expose the product to prohibited environmental conditions (see page 9).
- ▶ Do not use the product if it is damaged or in a questionable condition. If necessary, have it inspected, repaired or replaced by the manufacturer or a specialist workshop.
- ▶ The product is intended for use by one person only and may not be reused by other persons.

Signs of changes in or loss of functionality during use

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

3 Intended use

3.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

3.2 Operating conditions

The product is recommended for:

- Mobility grade 1: indoor walker
- Mobility grade 2: restricted outdoor walker
- Mobility grade 3: unrestricted outdoor walker
- Mobility grade 4: unrestricted outdoor walker with particularly high demands
- The maximum approved body weight is specified in the technical data (see page 13).

3.3 Environmental conditions

Storage and transport
Storage temperature: -20 °C to +60 °C, relative humidity: 20 % to 90 %, no mechanical vibrations or impacts
Allowable environmental conditions
Temperature for use: -10 °C to +45 °C

Allowable environmental conditions
Moisture: relative humidity: 20 % to 90 %, non-condensing
Prohibited environmental conditions
Chemicals/liquids: fresh water, salt water, perspiration, urine, acids, soapsuds, chlorine water
Solids: highly liquid-binding particles (e.g. talcum), dust, sand

3.4 Product combinations

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Limited combination options for Ottobock components

The prosthetic foot generates high moments in the ankle area. Use structural components with higher weight limits:

Body weight [kg]	up to 75	76 to 100	from 101
Structural component weight limit [kg]	≥ 100	≥ 125	≥ 150

3.5 Lifetime

Depending on the user's activity level, the maximum lifetime of the product is 3 years.

4 Scope of delivery

- Prosthetic foot, protective sock, instructions for use

5 Accessories

Designation	Reference number
Protective sock	SL=Spectra-Sock-7
Footshell with connection cover	2C5*
Connection cover	2C10*
Changing device	2C100

6 Preparing the product for use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- Observe the alignment and assembly instructions.

NOTICE

Altering the prosthetic foot or footshell

Premature wear resulting from damage to the product

- ▶ Do not alter the prosthetic foot or footshell.

INFORMATION

A plastic part protects the connecting section of the product from scratches during the alignment and fitting.

- ▶ Remove the plastic part before the patient leaves the fitting area.

6.1 Applying/removing the footshell

- ▶ Pull the protective sock over the prosthetic foot to prevent noises in the footshell.
- ▶ **CAUTION! Always use the prosthetic foot with the footshell.**
Apply or remove the footshell as described in the footshell instructions for use.

6.2 Bench Alignment

- > **Required materials:** goniometer 662M4, heel height measuring device 743S12, 50:50 gauge 743A80, alignment apparatus (e. g. PROS.A. Assembly 743A200)
- > **Transfemoral prosthesis or knee disarticulation prosthesis:**
Observe the instructions for use for the prosthetic knee joint.
- 1) Position the prosthetic foot in the alignment apparatus.
 - **Heel height: effective heel height of everyday shoe + 5 mm**
- 2) **Sagittal plane:** Align the prosthetic foot with the alignment reference line.
 - Position the alignment reference line in a posterior position in relation to the middle of the foot (**10 mm to 15 mm** for size 24 to 25, **15 mm** for size 26 to 28, **15 mm to 20 mm** for size 29 to 31)
- 3) **Frontal plane:** Align the prosthetic foot with the alignment reference line.
 - Set the exterior rotation to **approx. 5°** (alignment reference line in the middle of the second toe).
- 4) Install the structural components between the prosthetic foot and prosthetic socket.
- 5) **Sagittal plane:** Align the prosthetic socket with the alignment reference line.

- Determine the centre of the prosthetic socket proximally and distally with the 50:50 gauge and mark the centre line.
 - Mark the socket reference point on the centre line: **For transtibial prostheses, at the height of the centre of the patella.** The alignment reference line runs vertically through the socket reference point.
 - Adjust the socket flexion by rotating around the socket reference point: **individual residual limb flexion + 5°**
- 6) **Frontal plane:** Align the prosthetic socket with the alignment reference line.
- Take the abduction or adduction of the residual limb into account.

6.3 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the 3D L.A.S.A.R. Posture and adapting it as needed.

Optionally, the alignment recommendations (transfemoral modular leg prostheses: **646F219***, transtibial modular leg prostheses: **646F336***) can be requested from Ottobock as needed.

6.4 Dynamic Trial Fitting

The optimal gait pattern is developed during dynamic fitting. The prosthesis is aligned optimally in the frontal plane and the sagittal plane for this purpose.

- ▶ **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal planes is achieved after the leg begins to bear weight following the heel strike. Avoid medial movement of the knee joint.
 - If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction.
 - If the knee joint moves in the medial direction in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.
- ▶ After completing the dynamic fitting and the walking exercises: Remove the plastic part from the pyramid.

7 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

8 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.

- ▶ During regular inspection: Inspect the prosthesis for signs of wear and check its function.
- ▶ Conduct annual safety inspections.

9 Disposal

Do not dispose of the product in household waste. Improper disposal is harmful to health and the environment. Observe the instructions of the local authority regarding return, collection and disposal.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

11 Technical data

Size [cm]	24	25	26	27	28	29	30	31
Heel height [mm]	9 ± 5							
System height [mm]	18							
Product weight [g]	270	280	290	305	320	330	345	355
Max. body weight [kg]	Mobility grade 1 to 3: 136 kg, mobility grade 4: 100 kg							
Mobility grade	1 to 4							

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